

Your Life Science  
Compliance  
Partner

How to Improve Cleanroom Compliance

# A Proactive Checklist for Control & Predictability

**e|lab**

## Purpose

Cleanrooms are among the highest-stakes environments in the life sciences industry: even the smallest issue can lead to significant consequences. As regulatory expectations change and evolve, and as operational risks become increasingly complex, improving cleanroom compliance requires more than the traditional approach: It's time to adopt a more proactive, holistic, and data-driven strategy.

The objective of this checklist is to provide life science professionals with actionable insights on how to maintain cleanroom compliance, enhance process efficiency, and leverage calibration, temperature mapping, and real-time monitoring to support regulatory readiness.

## Scope

This checklist has been designed to ensure compliance with essential cleanroom requirements, including ISO 14644-4:2022, ISO 14644-3, EU GMP Annex 1, ICH Q9, and the ISPE Good Practices Guide.

By adhering to these standards, life science companies can not only maintain high levels of cleanliness and safety while fostering a productive environment that supports the integrity of their processes and products. Implementing these practices is crucial for successful operations in regulated industries.

# 01.

## Develop Risk - Based Temperature Mapping Strategies

Temperature mapping is more than a one-time activity or just a box to check during qualification and then forget. In regulated cleanroom environments, neglecting ongoing monitoring can lead to serious risks: Compromised product quality, delayed releases, or worse, threats to patient safety.

### Checkpoints

---

Temperature mapping protocols are developed using Quality Risk Management (QRM) principles, as outlined in the ISPE Good Practice Guide: Temperature Controlled Storage Systems (2nd Edition). This ensures a science-based approach to mapping, maintaining compliance, and safeguarding product integrity throughout the storage lifecycle.

---

Sensor placement is defined based on cleanroom size, equipment layout, HVAC dynamics, and process-critical zones.

---

Mapping is executed under both static and operational conditions to reflect true thermal behavior.

---

Heat maps and temporal trend analyses are used to identify spatial variability.

### Completed

---

## 01. Develop Risk - Based Temperature Mapping Strategies

### Checkpoints

Completed

Post-study mapping data is leveraged to refine sensor placement, EMS configuration, and risk mitigation strategies

Re-mapping intervals are defined based on environmental risk, seasonal variability, and room classification.

### PRO TIP

Use historical mapping data to validate EMS sensor placement, not just design intent.

# 02.

## Strategic Sensor Placement for Temperature Mapping

Temperature mapping is only as reliable as the accuracy - and justification - of your sensor placement. This section ensures that EMS sensor locations are chosen strategically and justified based on risk and regulatory expectations, resulting in high-resolution environmental awareness without redundancies or gaps.

### Checkpoints

### Completed

In small cleanrooms, a five-point sensor grid (four corners plus center) is deployed at a defined operational height.

In larger rooms, increased sensor density and vertical profiling are incorporated, with emphasis on high-activity and equipment-intensive zones.

Sensor locations are validated using historical mapping data to confirm environmental representativeness.

Mapping coverage is reviewed periodically to ensure alignment with the current facility layout and process evolution.

## 02.Strategic Sensor Placement for Temperature Mapping

### PRO TIP

Adjust your mapping strategy; even the duration of your study should adapt to the complexity of the room.



# 03.

## Improve Your Differential Pressure Mapping.

In GMP cleanrooms, the airflow matters just as much as temperature. Differential pressure is the difference in air pressure between two distinctive points in a system.

For a cleanroom environment, these two points are the inside of the cleanroom and its surrounding areas.

### Checkpoints

### Completed

Pressure differentials between cleanroom zones are mapped to confirm a minimum of 0.04" w.g. (10 Pa) as per EU GMP expectations.

Studies are conducted during dynamic operating conditions to capture real HVAC system behavior.

Variability, pressure reversals, and envelope leaks are identified and assessed against contamination risk.

Mapping results are used to adjust HVAC balancing, sensor setpoints, and alert thresholds.

### 03.Improve Your Differential Pressure Mapping.

**PRO TIP**

“Mapping is about confidence: Confidence that your space is performing as designed, not just when inspectors are watching, but every day, in every condition.”



# 04.

## Ensure Calibration Integrity Across Critical Instruments.

Calibration is a key factor in maintaining control in your cleanroom environments.

By repeatedly confirming your equipment performance, you ensure that all measurements remain accurate, traceable, and defensible under audit. A structured calibration program reduces uncertainty, strengthens data integrity, and supports timely decision-making in environments where precision is non-negotiable.

### Checkpoints

### Completed

All critical cleanroom instruments (e.g., differential pressure transmitters, particle counters, thermal probes) are included in a defined calibration schedule.

Calibration standards are traceable to national or international metrology references (e.g., NIST, DKD).

Calibration records are managed electronically with full audit trails, electronic signoffs, and change control.

Calibration tolerance limits are risk-assessed and aligned with process impact potential.

## 04.Ensure Calibration Integrity Across Critical Instruments

### Checkpoints

### Completed

Calibration tolerance limits are risk - assessed and aligned with process impact potential.

Calibration drift and out - of - tolerance trends are analyzed for early failure detection.

“If it wasn't documented, it did not happen!”

**Your calibration program should be well - designed and executed consistently:**

Consistency means you are following your approved SOPs.

Instruments and equipment should be classified according to their criticality (where the highest criticality receives the greatest scrutiny).

Calibration specifications, including test points, tolerances, and intervals, need to be well defined, applicable, and appropriate for the process they are controlling

# 05.

## An EMS Design for Daily Compliance.

Environmental monitoring systems (EMS) are only as effective as their design. To maintain compliance and prevent contamination, your EMS must reflect the true environmental behavior of the cleanroom, day in and day out.

This section outlines how to build an EMS that delivers clarity, accountability, and control.

### Checkpoints

### Completed

EMS sensor placement is guided by mapping outcomes, zone classification, and contamination risk levels.

A combination of in-room and in-duct sensors is implemented to ensure holistic environmental visibility.

Display units and dashboards are installed near cleanroom entries or control rooms for operator awareness.

Alarm logic is configured to prioritize actionable events while mitigating alert fatigue.

# 06.

## Advanced Monitoring Technologies Built for the Future.

As cleanroom operations evolve, supporting technology must follow. Monitoring platforms should offer flexibility, validated performance, and robust infrastructure, ensuring that data remains reliable, secure, and accessible, even in the face of change.

### Checkpoints

### Completed

A validated EMS platform following Annex 1.1 guidelines is deployed with features such as Power-over-Ethernet, real-time alerting, and ISO 14644-compliant reporting.

Sensor type (wired or wireless) is selected based on risk evaluation, cleaning compatibility, and facility infrastructure.

System resilience is ensured through backup data storage, two-way communication, and power-failure contingencies.

Periodic software updates and revalidations are performed under change control.

## 06. Advanced Monitoring Technologies Built for the Future.



### PRO TIP

Don't wait until validation to add or change monitoring locations: Going back and making facility changes may trigger the need to remap that space.

# 07.

## Maintain Regulatory Alignment & Data Integrity Across Systems.

Data collected from your EMS is only as valuable as its integrity. Clear access control, audit-ready documentation, and adherence to electronic record regulations provide the foundation for defensible, inspection-ready systems that support continuous compliance.

### Checkpoints

### Completed

Monitoring systems are implemented in compliance with ISO 14644, 21 CFR Part 11, EU Annex 11, and cGMP documentation principles.

Access control, data logging, time stamping, and audit trails are enforced.

All changes, overrides, and alarm acknowledgments are electronically tracked and reviewable.

Data integrity assessments are conducted routinely as part of internal QA or CSV review cycles.

# 08.

## Empower Your Cleanroom Staff.

Well-trained personnel are essential to cleanroom success. From following SOPs to responding to alerts, empowering your teams with knowledge, clarity, and cross-functional collaboration strengthens your quality culture and safeguards day-to-day control.

### Checkpoints

### Completed

Personnel working in the cleanroom know and follow SOPs.

Personnel go through recurrent training to ensure in-depth familiarity with the procedures and processes.

A cross-functional collaboration between QA, validation, maintenance, and operations is encouraged.

# 09.

## Use Trend Analysis To Guide Requalification & Continuous Improvement.

Environmental monitoring isn't just about capturing data: It's about how you use it.

Long-term trend reviews reveal patterns, flag anomalies, and guide timely actions, helping you stay in control of any issues that might arise and keep cleanroom performance within defined limits.

### Checkpoints

### Completed

EMS data is routinely compared to baseline mapping results to detect environmental drift.

Trend reviews are scheduled at defined intervals and after significant events (e.g., maintenance, layout changes, HVAC interventions).

Root cause investigations are triggered when environmental trends deviate from control expectations.

Re-mapping or requalification is initiated when thresholds of variability or non-conformance are met.

## 09. Use Trend Analysis To Guide Requalification & Continuous Improvement.

### PRO TIP

Once the data's in, don't just log and store it: Use it to drive long-term decisions like:

Where to place permanent EMS probes

---

Whether airflow needs to be balanced

---

Where unexpected risks are emerging over time

By completing this cleanroom compliance checklist and addressing these key focus areas, you're taking meaningful steps to maintain control and enhance predictability in your cleanroom and monitored spaces.

A holistic strategy combined with a proactive, data-driven approach is essential for building and sustaining a compliant environment. When you integrate risk-based mapping, routine calibration, intelligent monitoring, and empowered teams, you create a strong foundation for continuous compliance, reduced variability, and long-term operational excellence.